

PathScribe AI · ForMedrix

v0.9.0

# PathScribe by ForMedrix AI

## User Guide

Clinical Reporting System

For Pathologists and Residents

May 2026 · Confidential — PathScribe AI

---

|                                                |           |
|------------------------------------------------|-----------|
| <b>Table of Contents</b>                       | <b>2</b>  |
| <b>1. Getting Started</b>                      | <b>4</b>  |
| 1.1 Logging In                                 | 4         |
| 1.2 Biometric Setup                            | 4         |
| 1.3 The Navigation Bar                         | 4         |
| <b>2. The Home Screen</b>                      | <b>5</b>  |
| <b>3. The Worklist</b>                         | <b>6</b>  |
| 3.1 Status Tiles                               | 6         |
| 3.2 Table Columns                              | 6         |
| 3.3 Sorting                                    | 6         |
| 3.4 Pool Cases                                 | 6         |
| 3.5 Delegated Cases                            | 7         |
| <b>4. Opening a Case</b>                       | <b>8</b>  |
| 4.1 Navigation within a Case                   | 8         |
| 4.2 Patient History                            | 8         |
| <b>5. Diagnostic Information</b>               | <b>9</b>  |
| <b>6. Synoptic Reporting &amp; AI Copilot</b>  | <b>10</b> |
| 6.1 AI Copilot Suggestions                     | 10        |
| 6.2 Accepting and Rejecting Suggestions        | 10        |
| 6.3 Protocol Change Review (After Microscopic) | 10        |
| 6.4 Coding                                     | 10        |
| <b>7. Computational Sidecar</b>                | <b>12</b> |
| 7.1 Worklist Computational Icons               | 12        |
| 7.2 Synoptic Report — Computational Tab        | 12        |
| 7.3 Cancelling a Pending Order                 | 12        |
| 7.4 Ordering Additional Tests                  | 12        |
| <b>8. Report Preview</b>                       | <b>13</b> |
| <b>9. Case Finalisation</b>                    | <b>14</b> |
| 9.1 Pre-Finalisation Review                    | 14        |
| 9.2 After Finalisation                         | 14        |
| <b>10. Case Delegation</b>                     | <b>15</b> |

|                                                                      |           |
|----------------------------------------------------------------------|-----------|
| 10.1 Delegating a Case . . . . .                                     | 15        |
| 10.2 Receiving a Delegated Case . . . . .                            | 15        |
| <b>11. Case Search . . . . .</b>                                     | <b>16</b> |
| 11.1 Smart Identifier Box . . . . .                                  | 16        |
| 11.2 Search Filters . . . . .                                        | 16        |
| <b>12. Internal Messaging . . . . .</b>                              | <b>17</b> |
| 12.1 Reading Messages . . . . .                                      | 17        |
| 12.2 Composing . . . . .                                             | 17        |
| 12.3 Request Informal Review . . . . .                               | 17        |
| <b>13. Voice Commands . . . . .</b>                                  | <b>18</b> |
| <b>14. Your Profile . . . . .</b>                                    | <b>19</b> |
| <b>15. Quick Reference . . . . .</b>                                 | <b>20</b> |
| <b>Appendix A — Colour &amp; Visual Language Reference . . . . .</b> | <b>21</b> |
| A.1 Case Status Indicators (Worklist dots) . . . . .                 | 21        |
| A.2 Case Flags . . . . .                                             | 21        |
| A.3 TAT Progress Bars . . . . .                                      | 21        |
| A.4 Severity Badges (Quality Tab) . . . . .                          | 21        |

---

## 1.1 Logging In

Open Chrome or Edge and navigate to your PathScribe by ForMedrix AI address. Enter your email and password on the ForMedrix login screen and click Sign In. Credentials are issued by your system administrator.

## 1.2 Biometric Setup

On first login PathScribe may prompt you through the Biometric Setup Wizard (fingerprint or face recognition). You may skip this; it can be completed any time from your profile. Once enrolled, biometric sign-out replaces the password prompt in the Pre-Finalisation Review.

## 1.3 The Navigation Bar

| Element                    | Action                                                                    |
|----------------------------|---------------------------------------------------------------------------|
| PathScribe logo            | Click anywhere to return to the Home screen.                              |
| Enhancement Request (■)    | Submit feature requests or bug reports.                                   |
| Voice toggle (■)           | Enable/disable hands-free voice commands.                                 |
| Messages (■)               | Open the internal messaging drawer.                                       |
| Clinical Links (■)         | Quick access to CAP protocols, WHO Classification, and PathologyOutlines. |
| Your name / initials badge | Profile card — User Guide, Admin Guide (if admin), System Information.    |
| Sign Out (arrow)           | Log out of PathScribe.                                                    |

| Tile            | What it opens                                                        |
|-----------------|----------------------------------------------------------------------|
| Worklist        | Your active case queue — the primary daily view.                     |
| Contribution    | Personal productivity dashboard: workload, quality, and TAT metrics. |
| Search          | Cross-institutional case search with advanced filters.               |
| System Audit    | Audit trail of all system events (admin-visible).                    |
| My Contribution | Detailed personal metrics by specialty, TAT trend, and quality.      |
| Configuration   | System settings, user management, template builder (admin role).     |

The Worklist is your primary clinical view. STAT / critical cases are pinned at the top under a red Urgent banner; remaining cases follow in sort order.

### 3.1 Status Tiles

| Tile            | Colour | What it shows                                                                           |
|-----------------|--------|-----------------------------------------------------------------------------------------|
| Total Cases     | Grey   | Your full non-pool caseload. Always a summary — not a filter.                           |
| Pool Cases      | Indigo | Unassigned cases in the shared pool. Click a pool case to open the Pool Claim workflow. |
| Delegated to Me | Cyan   | Cases routed to you by another pathologist for review or second opinion.                |
| Critical        | Red    | STAT-priority or critically flagged cases.                                              |
| In Progress     | Blue   | Cases actively being worked on in the reporting workflow.                               |
| Needs Review    | Amber  | Cases with an AI-generated synoptic draft awaiting your sign-off.                       |
| Amended         | Purple | Finalised reports re-opened for amendment.                                              |
| Completed Today | Green  | Reports finalised today.                                                                |

### 3.2 Table Columns

| Column         | Contents                                                                                 |
|----------------|------------------------------------------------------------------------------------------|
| ■ (dot)        | Pulsing red dot = STAT/critical. No dot = routine.                                       |
| Case           | Accession number. Red = critical, cyan = routine. Click to open.                         |
| Patient        | Last, First name (PHI). Sortable alphabetically.                                         |
| MRN            | Medical Record Number (PHI-protected).                                                   |
| Sex            | M / F / O. Sortable.                                                                     |
| DOB (Age)      | MM/DD/YYYY with compact age label (e.g. 52y).                                            |
| Specimen(s)    | Up to 3 specimen label chips. Hover for full description.                                |
| Accession Date | Date specimen received in the laboratory. Sortable.                                      |
| Physician      | Requesting/submitting physician. Sortable.                                               |
| Flag(s)        | Coloured flag chips. Solid border = case-level. Dashed = specimen-level. Hover for name. |
| Status         | Coloured dot: Draft · In Progress · Finalizing · Finalized · Amended.                    |

### 3.3 Sorting

Click any column header to sort. Click again to reverse direction. PathScribe supports multi-column sort — up to three sort keys simultaneously. Sort preferences persist across sessions.

### 3.4 Pool Cases

Pool cases are unassigned cases visible to all pathologists at your institution. Click a pool case row to open the Pool Claim modal. Click Claim to assign it to yourself.

### 3.5 Delegated Cases

When another pathologist delegates a case to you, it appears in the Delegated to Me tile. Open the case and review the delegation note. You may accept the delegation or return the case with a comment.

---

Click any worklist row to open the case. The case layout has three panels:

| Panel  | Contents                                                                                               |
|--------|--------------------------------------------------------------------------------------------------------|
| Left   | Patient demographics, order details, clinical indication, requesting physician, and submitting client. |
| Centre | Diagnostic workspace: gross description, microscopic description, ancillary studies.                   |
| Right  | Synoptic reporting form, AI Copilot suggestions, and the Computational Sidecar tab.                    |

## 4.1 Navigation within a Case

Specimen tabs — multi-specimen cases show lettered tabs (A, B, C...) at the top of the right panel.

Breadcrumb navigation — click any crumb to return to that level. Previous / Next Case — use the arrow buttons in the case header bar (or voice commands).

## 4.2 Patient History

Click the Patient History link in the left panel to open the Patient History modal. This shows prior cases, relevant diagnoses, and clinical history for the patient. Dismiss with the close button, Esc, or click outside the modal.

The centre panel contains three autosaving rich-text fields. The AI Copilot reads these in real time and generates synoptic suggestions as you type.

| Field                   | Contents                                                                                                       |
|-------------------------|----------------------------------------------------------------------------------------------------------------|
| Gross Description       | Macroscopic findings — specimen dimensions, tumour measurements, margin distances, and macroscopic appearance. |
| Microscopic Description | Histological findings — tumour type, grade, invasion, margins, lymph nodes.                                    |
| Ancillary Studies       | IHC, FISH, molecular, flow cytometry, and any other supplementary results.                                     |

**Tip:** Enter the gross description first — the AI begins generating synoptic suggestions immediately. Complete or refine microscopic and ancillary text to trigger further suggestion updates.

PathScribe uses CAP-compliant and RCPATH synoptic templates. The template displayed matches the specimen type and selected protocol.

### 6.1 AI Copilot Suggestions

| Badge colour | Confidence | Meaning                                                                |
|--------------|------------|------------------------------------------------------------------------|
| Green        | ≥ 85%      | High confidence — strongly supported by diagnostic text.               |
| Amber        | 50–84%     | Moderate confidence — review the source citation carefully.            |
| Red          | < 50%      | Low confidence — treat as a starting point only; verify independently. |

### 6.2 Accepting and Rejecting Suggestions

| Action              | How                                                         |
|---------------------|-------------------------------------------------------------|
| Accept a suggestion | Click the tick (✓) icon — applied to the field immediately. |
| Reject a suggestion | Click the cross (✗) icon and enter your own value.          |
| Accept all          | Click Accept All at the top of the synoptic panel.          |

**Important:** AI suggestions are advisory only. You as the pathologist are responsible for reviewing and verifying every field before finalising. Confidence scores reflect computational certainty, not clinical accuracy.

### 6.3 Protocol Change Review (After Microscopic)

When you save a microscopic description, PathScribe AI re-evaluates whether the synoptic protocol assigned at gross is still appropriate. If the microscopic findings suggest a different protocol — for example, a core biopsy initially classified as Benign/NOS now shows invasive carcinoma — a Protocol Change Review modal appears.

**Important:** New in v0.9.0: The Protocol Change Review is raised automatically when AI detects a protocol mismatch after microscopic entry.

The modal shows each proposed change with:

- The current protocol and the proposed replacement side by side.
- The AI's clinical reasoning from the microscopic text.
- A confidence score for the proposal.
- A checkbox to accept or reject each individual change.

Select the changes you want to apply, or deselect any you want to keep as-is. Click **Apply Changes** to commit your selections. Click **Keep existing protocols** to dismiss without making changes. All decisions are recorded in the audit trail.

### 6.4 Coding

As synoptic fields are completed, PathScribe automatically suggests ICD-10 and SNOMED CT codes. Review and confirm, add manual codes, or override suggestions in the Coding panel.

The Computational Sidecar brings live laboratory results — IHC, molecular, FISH, cytogenetics — directly into your PathScribe workflow.

## 7.1 Worklist Computational Icons

| Icon colour | Status             | What it means                                                           |
|-------------|--------------------|-------------------------------------------------------------------------|
| Teal        | Pending            | Test ordered. No result received yet.                                   |
| Amber       | Preliminary        | Interim result available. Final result still pending.                   |
| Green       | Final              | Final result received. No immediate clinical action required.           |
| Red         | Final — Actionable | Final result received. Review required — clinical action may be needed. |

## 7.2 Synoptic Report — Computational Tab

Inside the Synoptic Report, click the Computational tab to open the full results panel. This shows specimen-grouped results with: result data table, status badge, specimen badge, AI concordance indicator, and AI extraction badge.

## 7.3 Cancelling a Pending Order

Click the trash icon beside a Pending test in the navigator. Read the confirmation dialog carefully — if the order originated in the LIS, you must also cancel it there.

## 7.4 Ordering Additional Tests

Click + Order additional test at the bottom of the Computational navigator. Stage tests for one or more specimens, then click Place N Orders to submit.

**Tip:** You can stage tests for multiple specimens before submitting — change the Apply To selector between batches and accumulate all orders before clicking Place Orders.

---

Click [Preview Report](#) in the case header to open the formatted print-ready report.

| Control            | Function                                                             |
|--------------------|----------------------------------------------------------------------|
| Page size selector | Choose A4, US Letter, US Legal, A3, or B5. Your preference is saved. |
| Margin controls    | Adjust top, right, bottom, left margins in mm. Saved per session.    |
| Ruler              | Shows margin boundaries above the document.                          |
| Print button       | Sends the report to your browser print dialogue.                     |

When all required synoptic fields are complete, click **Finalise**. Finalisation is a three-step gate:

| Step                        | What happens                                                                                                                                                                                                                    |
|-----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1 — Validation              | PathScribe checks all required fields. Missing fields are highlighted. The gate will not advance until all required fields are filled.                                                                                          |
| 2 — AI Triage Review        | Fields the AI flagged as low-confidence (below the configured threshold) are shown one at a time. Review and confirm, override, or skip each item. Keyboard shortcuts: Space/→ = Confirm, O = Override, S = Skip, Esc = Cancel. |
| 3 — Pre-Finalisation Review | Full-screen two-pane review before credentials are collected (see below).                                                                                                                                                       |

## 9.1 Pre-Finalisation Review

**Important:** New in v0.9.0: The Pre-Finalisation Review replaces the previous direct-to-password confirmation. It is mandatory for all sign-outs.

The Pre-Finalisation Review shows a full-screen two-pane layout before the pathologist signs:

- **Left pane — Specimen list:** All specimens and their synoptic reports in transmission order. Drag specimens to reorder them. Click Exclude on any synoptic to remove it from this transmission (it remains in PathScribe as a draft for a future addendum).
- **Right pane — Report preview:** Live Q&A; formatted preview of the complete report output as it will be transmitted to the LIS. Updates in real time as you reorder or exclude synoptics.
- **Signing panel:** Sign using biometric (fingerprint or face, if enrolled) or your login password. This constitutes your clinical authorisation for the transmission and is recorded in the audit trail.

## 9.2 After Finalisation

Once finalised, the case status changes to Finalized and the report is transmitted to the LIS (Co-Pilot mode) or locked in PathScribe (Standalone mode). The report is locked for editing. Amendments require re-opening by an administrator.

**Important:** Resident counter-sign workflow: If you are a Resident, your report is submitted for attending review rather than directly finalised. The attending must counter-sign to complete the process.

---

Delegation routes a case to a colleague for peer review, specialist opinion, second opinion, or subspecialty referral.

## 10.1 Delegating a Case

| Step | Action                                                                                                               |
|------|----------------------------------------------------------------------------------------------------------------------|
| 1    | Open the case and click the Delegate button in the case header.                                                      |
| 2    | Select the recipient from the staff directory (search by name or subspecialty).                                      |
| 3    | Choose the Delegation Type (Peer Review, Second Opinion, Subspecialty Referral, MDT Discussion).                     |
| 4    | Add a note explaining what you need from the recipient.                                                              |
| 5    | Click Confirm Delegation — the recipient receives a notification and the case appears in their Delegated to Me tile. |

## 10.2 Receiving a Delegated Case

Click the Delegated to Me tile. Open the case — the delegation note and requesting pathologist are shown at the top of the right panel. Complete your review and either return the case or accept full ownership.

Case Search (Home → Search tile) finds cases across your institution using any combination of filters. All active filters combine with AND logic.

## 11.1 Smart Identifier Box

| What you type             | Detection                                  | Badge colour |
|---------------------------|--------------------------------------------|--------------|
| S26-4200                  | Accession number (matches your lab format) | Purple       |
| Smith, Jane or Jane Smith | Patient name                               | Green        |
| 123456 (numeric only)     | MRN / Hospital ID                          | Cyan         |
| abc123 (mixed)            | All fields — searches simultaneously       | Amber        |

## 11.2 Search Filters

| Filter            | Notes                                                                             |
|-------------------|-----------------------------------------------------------------------------------|
| Accession Date    | Date window. Quick-pick: 7d, 30d, 90d, 1yr. Default is last 30 days.              |
| Status            | Grossed, Awaiting Micro, Finalizing, Completed. Multi-select.                     |
| Priority          | Routine or STAT.                                                                  |
| Flags             | Click Browse. Multi-select. Flags are additive (OR logic within the flag filter). |
| Synoptic Protocol | Filter by CAP or institutional protocol. Click Browse for the full library.       |
| Pathologist       | Search by surname or hospital. Multiple pathologists can be selected.             |
| Specimen Type     | Multiple specimens allowed.                                                       |
| Diagnosis         | Free-text terms. Multiple terms allowed. Press Enter or + to add.                 |
| SNOMED CT         | Four axes: Morphology, Body Structure, Procedure, Specimen.                       |
| ICD Codes         | Search across all active coding systems (ICD-10, ICD-11, ICD-O).                  |

**Tip:** Save a "Daily STAT" search with a 7-day window and STAT priority for instant morning triage.

---

PathScribe's internal messaging system is for secure lab communication. Click the speech bubble icon in the navigation bar to open the messaging drawer.

## 12.1 Reading Messages

Messages are listed in the left panel. Unread messages are bold with a blue dot. Urgent messages appear highlighted in red. Click any conversation to read the full thread.

## 12.2 Composing

Click the compose icon at the bottom of the message list. Type a recipient name in the To: field. Enter subject and body. Toggle Urgent if required. Click Send.

## 12.3 Request Informal Review

To request an informal review from a colleague without transferring case ownership, open the case and use the Request Review option in the case header. Select a colleague, choose a review type (Informal Review, Consultation, Second Opinion, Clinical Observation), add an optional note, and send. The colleague receives a message with a direct link to the case.

**Note:** Request Informal Review does not transfer case ownership. It sends a message only. Use Case Delegation (Section 10) to formally transfer a case.

PathScribe is fully voice-enabled — every action available via mouse or keyboard can be triggered by voice. Click the microphone icon to activate. The system is context-aware: the same phrase does different things on different screens.

| Area            | Example commands                                                                          |
|-----------------|-------------------------------------------------------------------------------------------|
| Navigation      | Open worklist · Open search · Open messages · Go home · Go back                           |
| Case navigation | Next case · Previous case · Open case [accession number]                                  |
| Synoptic        | Accept suggestion · Reject suggestion · Accept all suggestions · Next field · Save report |
| Finalisation    | Finalise report · Submit for review                                                       |
| Messaging       | Open messages · Compose message · Send message · Mark read                                |
| Search          | Search cases · Clear filters · Export results                                             |
| System          | Open configuration · Sign out · Open voice settings                                       |

**Best results:** Use Chrome or Edge. Ensure your voice profile (EN-US or EN-GB) is set correctly in your profile settings. Speak clearly after the activation tone. PathScribe recognises 139 registered actions and learns unrecognised phrases permanently.

---

Click your name or initials badge (top right) to open the profile panel.

| Option             | Details                                                                                                 |
|--------------------|---------------------------------------------------------------------------------------------------------|
| User Guide         | Opens this document in a new browser tab.                                                               |
| Admin Guide        | Opens the Admin Guide (visible to Admin and Pathologist-Admin roles only).                              |
| System Information | Application version, AI provider and model, API connectivity status, browser and OS details.            |
| Copy to Clipboard  | Generates a formatted support report. Name and ID are redacted — safe to share with PathScribe support. |

| Action                               | How                                                                                      |
|--------------------------------------|------------------------------------------------------------------------------------------|
| Return to Home                       | Click the PathScribe by ForMedrix AI logo (top left)                                     |
| Open Worklist                        | Home → Worklist tile · voice: "open worklist"                                            |
| Open a case                          | Click any worklist row                                                                   |
| Navigate specimens                   | Click specimen tabs A, B, C... in the right panel                                        |
| View patient history                 | Patient History link in the left panel                                                   |
| Accept AI suggestion                 | Click tick (✓) next to the suggestion                                                    |
| Accept all suggestions               | "Accept All" button at the top of the synoptic panel                                     |
| Review protocol change               | Protocol Change Review modal — appears automatically after microscopic is saved (v0.9.0) |
| View computational results           | Click Computational tab in the synoptic panel                                            |
| Order an additional test             | + Order additional test at the bottom of the Computational navigator                     |
| Preview report                       | Preview Report button in the case header                                                 |
| Finalise a case                      | Finalise button → 3-step gate → Pre-Finalisation Review → sign                           |
| Reorder synoptics before sign-out    | Drag specimens in the Pre-Finalisation Review left pane (v0.9.0)                         |
| Exclude a synoptic from transmission | Click Exclude in the Pre-Finalisation Review left pane (v0.9.0)                          |
| Delegate a case                      | Delegate button in the case header                                                       |
| Request informal review              | Request Review option in the case header                                                 |
| Claim a pool case                    | Pool Cases tile → click a pool case row → Claim                                          |
| Search all cases                     | Home → Search tile                                                                       |
| Send a message                       | Nav bar speech bubble → compose icon                                                     |
| Get help                             | Your initials badge (top right) → User Guide                                             |
| Support report                       | Your initials badge → System Information → Copy to Clipboard                             |

This appendix explains the meaning of every colour, icon shape, and status indicator used across PathScribe.

## A.1 Case Status Indicators (Worklist dots)

| Colour                          | Meaning                                                                                  |
|---------------------------------|------------------------------------------------------------------------------------------|
| Red — STAT / Urgent             | Case flagged STAT or has a pending urgent action. Requires immediate attention.          |
| Amber — Action Required         | Needs user input: outstanding query, delegation pending, or a required field is missing. |
| Blue — In Progress              | Actively being worked on (grossing, microscopy, or synoptic in draft).                   |
| Green — Complete / Signed Out   | Report finalised and countersigned where required.                                       |
| Grey — Awaiting Micro / Grossed | Specimen received and grossed; sections processing; no pathologist action yet required.  |
| Purple — Delegated              | Case delegated to another pathologist or in a review pool.                               |

## A.2 Case Flags

| Colour | Flag class          | Examples                                                             |
|--------|---------------------|----------------------------------------------------------------------|
| Teal   | Administrative      | RUSH, Frozen Section, Tumour Board, MDT Discussion, QI Review        |
| Cyan   | Computational       | IHC panels, molecular assays, BRAF, PD-L1, ER/PR/HER2                |
| Amber  | Quality / Review    | Amended Report, Addendum Added, Correlation Required, Second Opinion |
| Green  | Positive / Complete | AI coding applied, countersign received, research consent captured   |
| Red    | High-severity       | Critical findings or urgent process flags                            |

## A.3 TAT Progress Bars

| Colour | % of SLA consumed | Action                                                          |
|--------|-------------------|-----------------------------------------------------------------|
| Green  | < 70%             | Well within the SLA. No intervention needed.                    |
| Amber  | 70–90%            | Approaching the SLA limit. Monitor closely.                     |
| Red    | > 90%             | At or exceeding the SLA. Escalation or prioritisation required. |

## A.4 Severity Badges (Quality Tab)

| Badge       | Meaning                                                                             |
|-------------|-------------------------------------------------------------------------------------|
| LOW (green) | Minor documentation issue or concordant variance. No immediate escalation required. |

| Badge             | Meaning                                                                                                       |
|-------------------|---------------------------------------------------------------------------------------------------------------|
| MEDIUM<br>(amber) | Moderate discordance (e.g. grade upgrade) or process deviation. Discuss at QI review.                         |
| HIGH (red)        | Significant diagnostic discordance or patient-safety relevant amendment. Requires formal root-cause analysis. |